



ENDYNE, INC.

Laboratory Services

160 James Brown Drive
Williston, Vermont 05495
Ph 802-879-4333
Fx 802-879-7103

SAMPLING INSTRUCTIONS

Fats, Oils & Grease
Hexane Extractable Material

and/or

Total Petroleum Hydrocarbons
Silica Gel Treated - Hexane Extractable Material
SGT-HEM

Samples must be kept on ice until delivery at the laboratory

EPA Method 1664A

2--1 Liter Amber Bottles

Hold Time: 28 Days

If Both FOG and HEM-SGT is required, 4-1L Amber Bottles must be submitted

-Fill bottle to just below the neck and cap. Sample acidification is done with 3.0mls **1:1 HCl** to a pH < 2. Sample bottles are shipped with the acid inside. Caution should be used to prevent spilling the acid. Do not rinse bottle before filling.