



ENDYNE, INC.

Laboratory Services

Lebanon, New Hampshire
Plattsburgh, New York
Williston, Vermont
www.endynelabs.com

IOCs, NITRATE and/or NITRITE, RADIOLOGICAL

**CHLORINATED or NON-CHLORINATED
PUBLIC WATER SYSTEMS**

SAMPLING INSTRUCTIONS

Fill all bottles to the bottom of the neck of the Bottle

INORGANIC CONTAMINANTS (IOC)

Metals, including Uranium 1--16oz Plastic Hold Time: 6 Months
Acid Preservation at lab

Fluoride 1 --2oz Plastic Hold Time: 28 Days
Does not need Ice

NITRATE and/or NITRITE

Nitrate 1 --2 or 8oz Plastic Hold Time: **48 hours**
Must be Less than 6° Celsius

Nitrate & Nitrite 1 --2 or 8oz Plastic Hold Time: **48 hours**
Must be Less than 6° Celsius

RADIOLOGICAL

(Analysis performed by certified sub-contracted Laboratory)

Radium 226 &228 2 -- 1L Plastic Hold Time: 1 Month
Acid Preservation at lab
Does not need Ice

Gross Alpha 1 -- 16oz Plastic Hold Time: 1 Month
Acid Preservation at lab
Does not need Ice

Radon in Water 2 -- 40mL Glass Vials Hold Time **48hrs**
Must be delivered to Lab on Monday, Tuesday, or Wednesday **No Headspace**
Must be Less than 6° Celsius



